


# EURUS SERVICE KIT

## 1P07BJA0 EURUS WATER RELAY VALVE KIT




SPRING FOR PULL BACK VALVE

GUIDE FOR PULL BACK VALVE

DIAPHRAGM FOR PULL BACK VALVE

1S 3 1pce

1A0PJKA0  
DIAPHRAGM FOR AIR VALVE  
1S 3 1pce



GASKET M5  
1S 3 6pcs

CHECK VALVE M5 (BARB FITTING 1/8)  
for BIEN AIR MX2/MCX  
1S 3 2pcs

CHECK VALVE M5 (BARB FITTING 1.8)  
for Turbine/Scaler  
1S 3 4pcs

1A0PURA0  
SPONGE FOR OIL MIST  
1 1S 3 1pce

1E09J6A0  
PN8136 WASHER FOR WATER BOTTLE  
1 1S 3 1pce

1E09J4A0  
PN3072 VALVE ASSY FOR DCI SYRINGE  
1S 1set

1E09J5A0  
PN2230 O-RING  
1S 1set(12pcs)

1A025DB0  
DRAIN STRAINER  
1 1S 3 1pce



## 1P055ZA0 SERVICE KIT WATER/AIR MAIN MASTER VALVE

MAIN MASTER VALVE SPRING  
O-RING P-16  
POPPET ASSY FOR MAIN  
MASTER VALVE  
DIAPHRAGM FOR MAIN  
MASTER VALVE  
1 1pce(Water) 1S 3 2pcs

1E097DA0  
WATER FILTER ELEMENT  
(100µ)  
1 1S 3 1pce

1E07K9A0  
AIR FILTER ELEMENT 40  
1 1S 3 1pce

1A0GCEA0  
SOLID COLLECTOR FILTER-CONE SHAPE  
1 1S 3 2pcs

Each parts are included in below kit

1 1 year kit 1S 1 year severe kit 3 3 years kit