


EURUS SERVICE KIT

1P07BJ** EURUS WATER RELAY VALVE KIT



SPRING FOR PULL BACK VALVE

GUIDE FOR PULL BACK VALVE

DIAPHRAGM FOR PULL BACK VALVE

1S 1pce

1A0PJK**
DIAPHRAGM FOR AIR VALVE
1S 1pce

GASKET M5
1S 6pcs

CHECK VALVE M5 (BARB FITTING 2.2)
for BIEN AIR MX2/MCX
1S 2pcs

CHECK VALVE M5 (BARB FITTING 1.8)
for Turbine/Scaler
1S 4pcs

1A0PUR**
SPONGE FOR OIL MIST
1S 1pce

1E09J6**
PN8136 WASHER FOR WATER BOTTLE
1S 1pce

1E09J4**
PN3072 VALVE ASSY FOR DCI SYRINGE
1S 1set

1E09J5**
PN2230 O-RING
1S 1set(12pcs)

1A025D**
DRAIN STRAINER
1S 1pce

1G0AA7** SERVICE KIT WATER/AIR MAIN MASTER VALVE

MAIN MASTER VALVE SPRING

O-RING P-16

POPPET ASSY FOR MAIN
MASTER VALVE

DIAPHRAGM FOR MAIN
MASTER VALVE

1S 2pcs

1G0AZG**
WATER FILTER ELEMENT SET
(100μ)
1S 1pce

1G0AZH**
AIR FILTER ELEMENT SET
(5μ)
1S 1pce

1A0GCE**
SOLID COLLECTOR FILTER-CONE SHAPE
1S 2pcs

Each parts are included in below kit

1S 1 year severe kit