

Belmont Handpieces

Belmont Handpieces



<https://dental.takarabelmont.co.jp/>



Belmont Handpieces - your complete dental solution

With the Japan-made:

- Enhance your experiences at work
- Confirm your confidence in quality
- Improve reliability on the results



Symbol

Handpiece, air motor with illumination

This product is autoclavable up to max.135°C.

Coupling with LED illumination

This product can be washed via Thermo Disinfector.

Clean Head System

Air Turbines

Miniature Head

DA-390X Turbine

- Stainless Steel Body



DA-391X Turbine

- Stainless Steel Body



DA-312 Turbine

- Stainless Steel Body



DA-520 Turbine

- Stainless Steel Body



DA-385 LED Coupling

- Titanium Coating
- Stainless Steel Body
- Non-Retraction Valve



DA-525 Coupling

- Stainless Steel Body
- Non-Retraction Valve



DA-526 Coupling

- Non-Retraction Valve



Air Motors

DA-281 Air Motor

- Torque 1.5 N-cm and more



DA-290X-LC Contra-Angle

- 1:1 Direct Drive
- Stainless Steel Body



DA-290X-S Straight Handpiece

- 1:1 Direct Drive
- Stainless Steel Body



Category	Model	Type	Optic	Power	Spray	Speed (min-1)	Head Size	Bur Type	Chuck	Bearing	Weight (g)
Air Turbine	DA-390X	Miniature Head	LIGHT	23W	4x	380,000-450,000	Ø 10.6 x H12.4 mm	ø1.59-1.60mm Std	Push Button	Ceramic	Approx. 46
	DA-391X	Standard Head	LIGHT	26W	4x	325,000-430,000	Ø 12.1 x H13.3 mm	ø1.59-1.60mm Std	Push Button	Ceramic	Approx. 47
	DA-312	Standard Head	NO LIGHT	20W	4x	360,000-430,000	Ø 12.1 x H13.3 mm	ø1.59-1.60mm Std	Push Button	Ceramic	Approx. 45
	DA-520	Standard Head	NO LIGHT	17W	4x	350,000-430,000	Ø 12.1 x H13.3 mm	ø1.59-1.60mm Std	Push Button	Ceramic	Approx. 51
Air Motor	DA-281	Internal Water Spray	NO LIGHT	-	-	22,000±10%	-	-	-	-	Approx. 64
	DA-290X-LC	1:1 Contra-Angle	LIGHT	-	1x	Max 40,000	Ø 9.9 x H13.8 mm	Ø2.35 CA	Push Button	-	Approx. 78
	DA-290X-S	1:1 Straight Handpiece	NO LIGHT	-	1x	Max 40,000	-	Ø2.35 HP (CA*)	-	-	Approx. 79

* When using the bur stopper (Option).

JAPAN ● MADE